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Instruction Manual

DDpmmaCAST

DDpmmaCAST

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Material-related characteristics

Chem properties	100% (polymethylmetacrylate)
Density	1,19 g/cm ³
Color	ivory

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Explanation of symbols

Manufacturer	Consult instructions for use
Batch code	Catalogue No.
Store dry	Protect from direct sun light

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Instruction Manual

- General information**
DD pmma CAST discs can be machined with any milling system. Tungsten carbide tools for PMMA are required for machining.
- Storage**
Protect from heat and direct sunlight.
- Indication**
DD pmma CAST discs are not certified for intraoral use and therefore must not be used in the patients mouth. The material is manufactured exclusively for the dental casting technique.
- Finishing**
The milled frameworks can be separated from the blanks with suitable tungsten carbide cutters or cut-off wheels.
 - Spruing**
Please check the rules for your casting alloys or pressable ceramics to design sprues and runners for casting. Acrylic molds for casting should be covered with a thin layer wax before embedding.
 - Embedding**
Use a ring free embedding system or coat the ring with a ceramic ring liner before filling with investment. Mix the embedding material using vacuum. Follow the manufacturers instructions for use of the investment material, the alloy and / or press ceramic used.
- Disposal**
DD pmma CAST discs are made from pure PMMA. Waste can be recycled with other acrylics or can be disposed of through domestic residual waste system.