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Instruction Manual

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Material-related characteristics

Chemical composition [wt. %]

Polyethylenterephthalat-Glycol (PETG)	100
---------------------------------------	-----

Physical specifications

Flexural strength	[ISO 20795-1]	[MPa]	≥ 58
Flexural modulus	[ISO 20795-1]	[MPa]	≥ 1370
Fracture toughness	[ISO 20795-1]	[MPa m ^{1/2}]	≥ 3
Water absorption	[ISO 20795-1]	[µg/mm ³]	6
Solubility	[ISO 20795-1]	[µg/mm ³]	< 0.25

Explanation of symbols

Manufacturer	Date of manufacture	Use-by date
Batch Code	Catalog number	Medical Device
Keep dry	Keep away from sunlight	
Content (Quantity)	Consult instructions for use	

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

1. Indications

DD Bio Splint FLEX are transparent dental milling blanks made of PETG, for the manufacture of flexible splints, therapeutic splints, bite regulators and drilling templates, for short-term use in the oral cavity of up to 30 days.

2. Contraindications

- intolerance to the contained components
- crowns, bridges, model casting with clasps, bars, narrow sublingual bars and transversal connectors

3. Safety information

Please take note of the information in the latest version of the safety data sheet. Avoid the inhalation of milling dust during processing. Wear gloves, protective goggles and a face mask to avoid skin irritation.

4. Handling and storage

The medical device may only be processed by trained personnel. Only use clean, defect-free material from undamaged packaging and ensure that the blanks are stored at room temperature, in their original cartons, in a dry location and protected from light. Avoid all contact with fluids.

5. Milling and further processing

Milling must be performed using a milling strategy that is aligned with the material.

Only use one-edged milling tools in the milling machine or wet grinding procedures!

Further processing:

Carbide cutters with FSQ tooling and black ring cutters are well-suited to the rough cutting for titanium. The fine post-processing can be performed using simple, cross-toothed HM cutters. The retaining bars are removed carefully and without high pressure. It is possible to polish the material; however, the contact pressure should be kept as low as possible since the resulting heat could otherwise result in smearing or melting.

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6. Cleaning

- ultrasonic bath (max. 40 °C) with a gentle cleaning agent
- Disinfection with IPA or med. Disinfectants
- Gamma and ETO sterilization possible

Do not use the following methods for cleaning:

- steam blasting
- autoclave method (e.g. sterilization / disinfection)
- organic solvents and acids (e.g. acetone, hydrochloric acid)
- highly concentrated alcohols
- chemical denture cleaners
- abrasive cleaning agents

7. Disposal

Remains of PETG blanks can be disposed of in the regular household garbage.

Our devices are subject to strict quality controls and reflect the current state of technology. We are continuously enhancing our devices and therefore reserve the right to make changes. The latest version of the instructions for use can be found on our website at: www.dentaldirekt.de/en.

This version replaces all previous versions.