


TECHNICAL DATA

Nacera® Pearl Q³ Multi-Shade

Type/ Color: A light, A dark, B light, B dark, C light

PHYSICAL PROPERTIES & TEST METHODS

| NORM | Property | Requirement | Typical Values |
|------|--------------------------|-----------------------------|--------------------------------------|
| | Optical properties | | Ultra high translucent |
| 6872 | Biaxial bending strength | ≥ 500 [MPa] | 600 |
| 6872 | Radioactivity | ≤ 1 [Bq/g] | < 0,0051 |
| 6872 | CTE | Determination | 9,8*10 ⁻⁶ K ⁻¹ |
| 6872 | Chem. Solubility | ≤ 100 [µg/cm ²] | < 13 |

Notified body: mdc GmbH  0483

These properties were measured on test samples. The values are typical material properties and may vary according to product configuration, geometry and manufacturing process.

Rev. 01-2021